

# Clinical research with the rtx1™ Adaptive Optics retinal camera

## Summary of published results in **age-related macular degeneration**

Dry AMD is the most common type of AMD, accounting for 90% of diagnosed cases. In this form of the disease, the breakdown of macular cells results in growing areas of geographic atrophy (GA). However, the progression is very slow and its detection usually takes months.

Thanks to Adaptive Optics (AO) technology, the rtx1™ enabled visualizing retinal changes in AMD patients at the cellular level. Moreover, since the rtx1 images are distortion-free, follow-up images could be aligned with high precision to study the dynamics of such changes over shorter time scales.

Clinical studies using the rtx1 have resulted in several new findings:


- In the early stages of AMD, rtx1 images showed that drusen cause changes in the contrast of the cone mosaic. Conventional drusen and reticular pseudodrusen were characterized by different reflectivity profiles<sup>1-3</sup>.
- On eyes with dry AMD, GA borders and spared foveal areas appeared more detailed on rtx1 images than on SLO and autofluorescence (AF) images<sup>4,5</sup>. Time-lapse rtx1 imaging enabled tracking displacements of GA borders with micrometer precision, and progression could be detected in less than 1 month<sup>4</sup>.
- The rtx1 also revealed a new candidate biomarker for dry AMD: the hyporeflective clumps (HRCs) which accumulate and migrate during disease progression. Although HRCs are invisible with other imaging techniques, the rtx1 enabled to observe their motion within a few days<sup>1,4</sup>.
- In addition to dry AMD, the rtx1 is also being used in two clinical investigations of stem-cell therapies for the wet form of late-AMD. rtx1 images helped assessing the survival of the implanted RPE cells<sup>6</sup> in one investigation, and of cone cells<sup>7</sup> in the other.

“The current technological level of robustness and the possibility to obtain quantitative biomarkers already permits the integration of AO in large scale trials in AMD

Paques et al., *Progress in Retina and Eye Research*, 2018

“Using AO NIR imaging, we could measure the progression of emerging atrophic areas at a relatively small temporal and spatial scale.

Gocho et al., *Investigative Ophthalmology & Visual Science*, 2013




Geographic atrophy progression monitored with the rtx1. Top: a nascent GA (arrowhead) is visibly growing. Bottom: enlarged part of the baseline image with the positions of the main GA border overlaid. Credit: Quinze-Vingts National Eye Hospital, Paris

# Clinical research with the rtx1™ AO camera

## Summary of published results in **age-related macular degeneration**

” Preservation of functional cone photoreceptors could be demonstrated on en face AO images in areas of foveal sparing that highlights the utility of this imaging modality in the evaluation of emerging treatments for GA.

Querques et al., *Retina*, 2016



Multimodal imaging of foveal sparing in a dry AMD case. Compared with the AF image, the rtx1 image shows the spared area (\*) more sharply and reveals HRCs in the atrophic area (dark dots). Credit: Quinze-Vingts National Eye Hospital, Paris


” Migration of HRCs is a highly dynamic process in AMD; it can indeed be detected over a timescale of days while atrophy progression is only detectable over a timescale of weeks.

Paques et al., *Progress in Retina and Eye Research*, 2018



www.imagine-eyes.com

18 rue Charles de Gaulle  
91400 Orsay, FRANCE  
+33 (0) 1 64 86 15 66  
contact@imagine-eyes.com



rtx1 images of conventional drusen (left) with hyperreflective rings, and of pseudodrusen (right) with hyporeflexive rings. Credit: Quinze-Vingts National Eye Hospital, Paris

” The rtx1 allowed us to directly observe stem cell-derived RPE cells after their transplantation in a patient's retina. Thanks to the rtx1's microscopic resolution, we could verify that the mosaic arrangement of these cells was similar to that of natural RPE cells, and stable over time.

Dr. Seiji Takagi, Kobe City Eye Hospital, Japan, 2019

### References

1. Paques, M. et al. Adaptive Optics Ophthalmoscopy: Application to Age-Related Macular Degeneration and Vascular Diseases. *Progress in Retinal and Eye Research* 66, 1-16 (2018)
2. Querques, G. et al. Appearance of Medium-Large Drusen and Reticular Pseudodrusen on Adaptive Optics in Age-Related Macular Degeneration. *British Journal of Ophthalmology* 98, 1522–1527 (2014)
3. Mrejen, S., Sato, T., Curcio, C. A. & Spaide, R. F. Assessing the Cone Photoreceptor Mosaic in Eyes with Pseudodrusen and Soft Drusen In Vivo Using Adaptive Optics Imaging. *Ophthalmology* 121, 545–551 (2014)
4. Gocho, K. et al. Adaptive Optics Imaging of Geographic Atrophy. *Investigative Ophthalmology & Visual Science* 54, 3673–3680 (2013)
5. Querques, G. et al. Adaptive Optics Imaging of Foveal Sparing in Geographic Atrophy Secondary to Age-Related Macular Degeneration. *Retina* 36, 247–254 (2016)
6. Takagi, S. et al. Evaluation of Transplanted Autologous Induced Pluripotent Stem Cell-Derived Retinal Pigment Epithelium in Exudative Age-Related Macular Degeneration. *Ophthalmology Retina*, doi:10.1016/j.oret.2019.04.021 (2019)
7. da Cruz, L. et al. Phase 1 Clinical Study of an Embryonic Stem Cell-Derived Retinal Pigment Epithelium Patch in Age-Related Macular Degeneration. *Nature Biotechnology* 36, 328–337 (2018)

(\*) Quantifications based on adaptive optics retinal images have not been approved for clinical use.

The rtx1 is an approved medical device in the European Union (device class 2a) and in Japan. In the USA, the rtx1 has not received FDA clearance. It is an investigational device and requires Institutional Review Board (IRB) oversight for use in any research application. Further information is provided in the user's documentation.