

# vitra 689™

NEXT  
GENERATION  
PDT LASER



Red 689 nm Laser  
for Photodynamic Therapy

# vitra 689™

## TECHNICAL SPECIFICATIONS

### VITRA 689

Laser Source:	Semiconductor laser diode
Wavelength:	689 nm +/- 3 nm
Output Power:	Variable to meet 600 mW/cm <sup>2</sup> at selected spot size
Maximum Power:	400 mW
Exposure Time:	83 s
Aiming beam:	650 nm
Size:	18 (H) x 19.5 (W) x 30 (D) cm 7.1" (H) x 7.7"(W) x 11.8"(D)
Weight:	5.6 kg – 12.3 lbs
Cooling:	Air
Power requirements:	100 to 240 VAC, 350 VA, 50/60 Hz


### SLIT LAMP ADAPTERS

Slit lamp compatibility:	Quantel Medical slit lamp, Haag Streit 900 BM & BQ and clones
Zoom:	Motorized
Spot size:	Continuously adjustable from 1200 to 5400 µm (Up to 8000 µm with a magnifying contact lens)
Weight:	0.9 kg – 1.98 lbs

Designed and Manufactured by Quantel Medical in France  
© 2019. Quantel Medical, Vitra 689® is trademark of Quantel Medical. All rights reserved.



[www.quantel-medical.com](http://www.quantel-medical.com)



0459



ISO 9001 : 2015 – ISO 13485 : 2016

**Headquarters**  
Quantel Medical  
11, rue du Bois Joli - CS40015  
63808 Cournon d'Auvergne – FRANCE  
Tel: +33 (0)4 73 745 745  
Email: [contact@quantel-medical.fr](mailto:contact@quantel-medical.fr)

**North America**  
Quantel Medical USA  
2313 Ridge Road Suite 105A  
Rockwall, TX 75087 – USA  
Tel: +1 888 660 6726  
Email: [info@quantelmedical.com](mailto:info@quantelmedical.com)

**Poland**  
Quantel Medical  
ul. Raclawicka 93  
02-634 Warszawa - POLAND  
Tel: +48 22 5210111  
Email : [info@quantel-medical.pl](mailto:info@quantel-medical.pl)

**Representative Offices**  
Thailand, Chiang Mai  
Brazil, Rio De Janeiro